Waist Massager
W9Pro

Please start and experience it now

Front of belt

Wear this side outwards

Back of belt

Wear this side inwards and adjust based on your waist

Message belt

On/Off button, press and hold for 1.5S to turn on/off

Type-C charging port

Infrared irradiation lamp

Indicator lights (battery level indicator lights)

When fully charged, all 4 lights are on, and the number of white indicator lights represents the remaining power.
Only one indicator light remaining on in red means low battery, and it is recommended to charge it in time.
(The wireframe is for illustration purposes only and may differ in appearance from the actual product.)

Operating Instructions

Step 1: After receiving the product, download the SKG Health APP (refer to the smart control procedures on APP). If you have purchased the remote control, you can also use the key to control the belt (please refer to the package insert of the remote control for how to use);
Step 2: Press and hold the On/Off button for 1.5 seconds to turn it on. After the electrode pads are warm, wear the massage belt around the waist, and then tighten the belt, so that the massage head electrode pads can be attached tightly on the skin;
Step 3: Open the mobile APP or use the remote control to control the intensity, mode and hot compress intensity level;
Step 4: After massage, the product automatically turns off. Remove the belt and store in a clean and safe place.

Product basic parameters

Product name	Waist Massager		
Model	SZ133YA	Rated charging voltage	5V =
Battery capacity	2450mAh/9.07Wh	Rated input power	10W

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SGK Health Technologies Co., Ltd.