



European
Funds
Smart Growth



Republic
of Poland

European Union
European Regional
Development Fund



RACER®

Special needs stroller

RACER® is a special needs stroller designed for children with limited mobility, providing stability, security, and proper seating support. It allows for both an upright and reclining position, making it ideal for daily use at home or in rehabilitation centers. The stroller helps parents with care, feeding, and outdoor activities, enabling children to spend time actively with their families. Its adjustable features ensure comfort, while the collapsible frame makes transport easy. The breathable, washable upholstery enhances durability and hygiene.



Light-
weight

24⁺

24-month
warranty



Stable
construction



Smooth and stepless
adjustment

Colors

Tapicerka ATMO™

Breathable and comfortable

HYDRO™ upholstery

Waterproof and durable



Grass



Sea



Lime



Aqua



Denim



Steel



Cobalt



Wheels with off-road tread guarantee the highest level of comfort even in difficult terrain.



Two-piece upholstery – much more comfortable and durable. Simple disassembly makes it easy to clean and maintain hygiene.



120 - 180 cm

Indicative user
height



35 - 90 kg

Max. user
weight



For more details about
the product, please visit
our website.



Innovative medical products with
a smooth and stepless adjustment

This is a medical device.
Use it in accordance with the user manual or the label.

INDICATIONS

Cerebral palsy (CP) – various types and forms | Brain and spinal cord injuries of various etiologies with limb paresis, severe balance disturbances, and coordination impairments | Demyelinating diseases (multiple sclerosis – MS, amyotrophic lateral sclerosis – ALS, and others) | Meningocele | Muscle diseases with paresis and paralysis (including dystrophy, myasthenia gravis, polymyositis) | Genetic syndromes with limb paresis | Degenerative and metabolic diseases of the nervous system | Degenerative and rheumatic joint diseases | Other diseases with limb paresis, paralysis, and musculoskeletal impairments | After injuries, following traumatic brain injuries or spinal cord injuries with limb paresis and paralysis, particularly with quadriplegia (as part of the rehabilitation process or as permanent orthopedic support)

CONTRAINDICATIONS

Observation: Monitor user during standing. Stop therapy if paleness, sweating, fainting, joint overload, or body deformations occur; consult with a doctor or physiotherapist | **Absolute Contraindication:** Critical condition of the user. | **Relative Contraindications:** Resting hypertension > 200/110 mmHg, pulmonary hypertension, heart rhythm disorders, unstable asthma, unstable metabolic issues, moderate heart valve disease. | Conditions like healing fractures or severe osteoporosis that prevent weight-bearing | Incomplete fracture healing (unless approved by a doctor) | Postoperative joint load restrictions | Pain from dislocation, excessive tissue stretch, or pressure; dislocation without pain is not contraindicated | Significant hip/knee flexion or ankle contractures; ensure no excessive pressure or stretch on affected areas.

Manufacturer: Akces-MED Sp. z o.o. | Poland, 36-002 Jasionka | Jasionka 955B